



Cost-Effective Interoperability

Build an unlimited number of flexible and seamless interfaces across healthcare systems

Scalable interoperability that streamlines the exchange of health information enables more efficient operations and promotes better outcomes across communities.

Secure, flexible, vendor-agnostic data exchange

Mirth® Connect by NextGen Healthcare allows you to optimize your interoperability needs regardless of the system, message type, and transfer protocol—with a low cost of ownership and a large user community. This proven integration engine has a template-driven architecture, providing you quick, easy, and flexible channel development. Mirth Connect is also a Java-based application enabling you to use your server preference for deployment. We provide the flexibility and resources your organization needs to elevate your performance, with the option of expert-level support, commercial extensions, and tiered pricing.

Connect to any system over any protocol

- FHIR
- HTTPS
- Web Services (SOAP)
- REST API
- MLLP
- TCP/IP
- Database
- File System
- Email
- FTP / FTPS / SFTP
- JMS
- DICOM
- Serial
- ASTM E1381

*Custom Connectors can also be written using Java or JavaScript

Flexibility with expertise

Founded as an open-source product, Mirth Connect is available via both open-source and through commercial licensing. Confidently develop, launch, and maintain your interface engine—backed by our expert support—and choose the level that works for you.

Mirth Connect—Support suited to your needs

	Silver <small>(Direct Healthcare Providers Only)</small>	Gold	Platinum
SSL Manager	✓	✓	✓
Channel history plug-in	✓	✓	✓
Message generator plug-in	✓	✓	✓
Multifactor authentication (DUO, TOTP)		✓	✓
Role-based access control		✓	✓
Email reader		✓	✓
Serial connector		✓	✓
FHIR R4		✓	✓
CDAPI java library		✓	✓
Interoperability connector suite		✓	✓
LDAP authentication		✓	✓
ASTM E1394 data type		✓	✓
ASTM E1381 transmission mode		✓	✓
Advanced clustering			✓
Advanced alerting			✓
24x7 emergency production support		✓	✓
Incidents	10	40	Unlimited

***OEM and Perpetual License Models available upon request**

NextGen Healthcare also offers commercially supported platforms to deploy Mirth Connect and an in-house professional services team for custom channel development support. If you're looking to license Mirth Connect, we offer a comprehensive OEM license that gives you all the tools you need.

“I’ve been using interface engines for about 25 years. Mirth Connect, by far, blows away the other three engines that I’ve worked with.”

Robert Werner
Senior Interface Developer
 University of Toledo Medical Center

BELIEVE IN BETTER.®

Visit nextgen.com/connect to request a demo.