



# Cost-Effective Interoperability

Build an unlimited number of flexible and seamless interfaces across healthcare systems

Scalable interoperability that streamlines the exchange of health information enables more efficient operations and promotes better outcomes across communities.

## Secure, flexible, vendor-agnostic data exchange

Mirth® Connect by NextGen Healthcare allows you to optimize your interoperability needs regardless of the system, message type, and transfer protocol—with a low cost of ownership and a large user community. This proven integration engine has a template-driven architecture, providing you quick, easy, and flexible channel development. Mirth Connect is also a Java-based application enabling you to use your server preference for deployment. We provide the flexibility and resources your organization needs to elevate your performance, with the option of expert-level support, commercial extensions, and tiered pricing.

## Connect to any system over any protocol

- FHIR
- HTTPS
- Web Services (SOAP)
- REST API
- MLLP
- TCP/IP
- Database
- File System
- Email
- FTP / FTPS / SFTP
- JMS
- DICOM
- Serial
- ASTM E1381

\*Custom Connectors can also be written using Java or JavaScript

## Gain more flexibility with expert support

Mirth Connect Premium goes beyond the open-source solution to help you develop, launch, and maintain your interface engine with our expert support.

## Mirth Connect Premium—Support suited to your needs

	Silver <small>(Only for NextGen® Enterprise clients)</small>	Gold	Platinum
SSL Manager	✓	✓	✓
Channel history	✓	✓	✓
Message generator	✓	✓	✓
Multifactor authentication (DUO, TOTP)		✓	✓
Role-based access control		✓	✓
Email reader		✓	✓
Serial connector		✓	✓
FHIR R4		✓	✓
CDAPI Java library		✓	✓
Interoperability connector suite		✓	✓
LDAP authentication		✓	✓
ASTM E1394 data type		✓	✓
ASTM E1381 transmission mode		✓	✓
Advanced clustering			✓
Advanced alerting			✓
24x7 emergency production support		✓	✓
Support cases	10	40	Unlimited

**\*OEM and Perpetual License Models available upon request. Mirth® Connect Cures Extension can be included with any Gold or Platinum Bundle for Mirth Connect versions 4.2 or higher.**

NextGen Healthcare also offers commercially supported platforms to deploy Mirth Connect and an in-house professional services team for custom channel development support. If you're looking to license Mirth Connect, we offer a comprehensive OEM license that gives you all the tools you need.

“I’ve been using interface engines for about 25 years. Mirth Connect, by far, blows away the other three engines that I’ve worked with.”

**Robert Werner**  
Senior Interface Developer  
University of Toledo Medical Center

# BELIEVE IN BETTER.®

Visit [nextgen.com/connect](https://nextgen.com/connect) to request a demo.